

Smart walker for clinical rehabilitation

Table of Contents

1. Problem Statement
2. Proposed Solution
3. How will we achieve the Specifications ?
4. Technical Details
5. Gantt Chart
6. Key Risks and Mitigate Strategies
7. Deliverables

01

Problem Statement



What is the problem?

Patients undergoing rehabilitation often use walkers incorrectly, leading to:

- Improper foot placement.
- Uneven weight distribution.

These issue causes :

- Increased risk of injuries.
- Slower recovery times.


Why is this Important ?

- Improper walker usage can lead to long-term mobility issues.
- There is a need for a smart, assistive device to improve rehabilitation outcomes

Current solutions lack real-time feedback to guide patients effectively !

02

Proposed Solution



Our Solution

A Smart Walker that provides real-time feedback to users during rehabilitation.

Key features:

- **Correct foot placement guidance** using a stereo camera.
- **Weight distribution monitoring** using load cells on each leg.
- **Real-time feedback** via haptic sensors, LCD display, and voice feedback.

Value Proposition:

- Reduces risk of injuries during rehabilitation.
- Accelerates recovery by ensuring proper walker usage.
- Provides clinicians with data to monitor patient progress.

03

How Will We Achieve the
Specifications?

Methodology and Component Selection:

- **Research & Literature Review:** Identified gaps in existing walker technologies.
- **Professor Meet:** Discussions with Prof. Neeta (BSBE Dept) to identify problems to analyse.
- **Cost Optimization:** Selected affordable components considering accuracy factor.



Technical Details


Sr. No	Usage (Name of the "Circuit block/Functional block" in which components is used)	Part description (OpAMP,ADC,DAC etc)	Manufacturing Part Number	Package type/Footprint	Available in WEL?	Vendor	Required Qty + Spare	Cost per each item	Total Cost
1	Weight Sensing Block	Load Sensor	613 N	load cell sensor	No	Robu	4+2	108	648
2	Weight Sensing Block	ADC	HX711	Breakout Board	No	Robu	4+2	40	240
4	Compute Block	Raspberry Pi 5	Raspberry Pi 5 Model 8GB	Development Board	No	Robu	1+0	8209	8209
5	Display Module	HDMI Capactive Touch	N/A	N/A	Yes		1+0		
6	Stereo Camera Module	Stereo Camera Module	IMX219-83	Camera Module	No	Robu	1+0	4799	4799
7	Haptic Feedback	DC Vibration Motor Module	N/A	Vibration Module	No	Robu	2+2	54	216
8	Audio System	Digital Amplifier Board	PAM8403	Digital Amplifier Board	No	Robu	1+1	65	130
9	Audio System	Speaker	N/A	8OHM Trumpet	No	Robu	2+2	42	168
10	Power Source	Power Bank	APB-15	Power Bank	No	Amazon	1+0	1899	1899
11	LED Feedback	LED	CJMCU-123 WS2812	LED Module	No	Robu	4+2	15	90
12	Data Logging	SD Card	Raspberry Pi Micro SD	Memory Card microSD	No	Robu	1+0	399	399
13	Walker	Walker			No	Amazon	1+0	853	853
								TOTAL	17651

Justification of Component Choices:

- **Load Sensor (613 N)**: Chosen for high accuracy and durability in measuring weight distribution on the walker legs.
- **ADC (HX711)**: Selected for its compatibility with load cells and ability to provide precise digital output for weight measurement.
- **Raspberry Pi 5 (8GB)**: It should provide sufficient processing power for real-time data processing and integration of multiple sensors.
- **HDMI Capacitive Touch Display**: Offers an intuitive user interface for real-time feedback and guidance.
- **Stereo Camera Module (IMX219-83)**: Chosen for accurate depth perception and cost-effectiveness compared to depth cameras.

Justification of Component Choices:


- **DC Vibration Motor Module:** Provides tactile feedback to users, enhancing the real-time guidance system.
- **Digital Amplifier Board (PAM8403):** Ensures clear audio feedback with low power consumption.
- **Speaker (8OHM Trumpet):** Delivers loud and clear voice feedback for user guidance.
- **Power Bank (APB-15):** Provides portable and enough power for the entire system.
- **LED Module (WS2812):** Used for visual feedback (e.g., color-coded alerts) to guide users.
- **SD Card (Raspberry Pi Micro SD):** Enables data logging for monitoring user progress over time.
- **Walker:** Acts as the physical base for integrating all components.




05

GANTT CHART

GANTT CHART





Key Risks and Mitigation Strategies

Risks:

- Load cells may provide noisy data.
- Integration challenges between hardware and software.
- Stereo Camera calibration error.
- Delay in Feedback.
- R-pi processing limitations.

Mitigations:

- Regularly calibrate load cells to maintain accuracy.
- Test each block individually before integration.
- Perform regular calibration to ensure accurate depth measurements.
- Optimize code for efficiency.

07

Deliverables

Expected Outcomes:

- A functional smart walker prototype.
- Real-time feedback system for foot placement and weight distribution.
- Data logging for clinician analysis.
- Documentation (technical report, user manual).

Impact:

- Improved rehabilitation outcomes for patients.
- Potential for commercialization in the healthcare industry.